

A Plan for the Development of an Integrated Common Diagnostics System

George Ludwig, Ph.D.
Diagnostic Systems Division
USAMRIID - Ft. Detrick

Program Framework

- Problem
- Requirement
- Solution



USAMRIID

THE PROBLEM

(Medical Perspective)

Diagnostics and Patient Care: Goals

- Reduce/prevent mortality
- Reduce/prevent morbidity
- Maintain unit readiness



USAMRIID


Impact of Diagnostics on Patient Care

- Immediate postexposure (up to 24 h)
 - ✓ very low concentration of agent
 - ✓ **IMPACT**
- Acutely ill
 - ✓ low concentration of agent
 - ✓ **IMPACT**
- Critically ill
 - ✓ High concentration
 - ✓ **IMPACT**



USAMRIID

Operational Scenario



USAMRIID

Diagnostic Selection Criteria

- Speed



USAMRIID

Classical Methods for Identifying Biological Agents

- Culture Isolation (1 - 30 days)
- Animal inoculation (2 - 30 days)
- Antigen detection (4 - 18 hours)
- Antibody detection (2 hr - 10 days)

Shipping - 24 -72 hrs!



USAMRIID

Diagnostic Selection Criteria

- Speed
- Sensitivity



USAMRIID

Classical Methods for Identifying Biological Agents

- Culture Isolation
 - High
- Animal inoculation
 - High
- Antigen detection
 - Low to Moderate
- Antibody detection
 - High




Diagnostic Sensitivity Requirements

Agent	Infective Dose	Agent	Infective Dose
Anthrax	8,000 to 50,000 spores	Smallpox	10-100 organisms*
Brucellosis	10-100 organisms	VEE	10-100 organisms
Plague	100-500 organisms	Viral Hemorrhagic Fevers	1-10 organisms
Q-fever	1-10 organisms	Botulinum Toxins	~70 ng
Tularemia	10-50 organisms	Staph Enterotoxin B	~30 ng




Anthrax in Swabs

Average inhaled dose = 1.7×10^5 Spores
(between 3 and 20 LD₅₀)




Operational Scenario



USAMRIID

Operational Scenario



Time!



USAMRIID


Diagnostic Markers

- Agent
- Antigen
- Antibody
- Nucleic Acid
- Cytokines/Chemokines
- Hormones
- Host Gene Expression




USAMRIID

Diagnostics Time Curve -- VEE



Symptoms



The Requirement

(DoD Perspective)

- Simultaneous identification ≥ 8 agents
- Identify biological agents at relevant conc.
- Sensitivity ≥ 0.98
- Specificity ≥ 0.98
- Approved or exempted by the FDA
- Upgradeable – equipment/assays
- Sample processing ≤ 20 min
- Total assay time ≤ 25 min
- Self-calibrating with failure alert



USAMRIID

The Requirement

(DoD Perspective)

- Protect and preserve samples
- Software/Communications compatibility
- Set up \leq 30 minutes
- Visual and audible positive identification signal
- Operable using all army power systems
- Battery operated \geq 12 hours
- Back-up power capable
- Man portable
- Day/Night Operations.
- Onboard data storage.




USAMRIID

Comprehensive Integrated Diagnostic System




Automated Sample Processing



Immunodiagnosics




Nucleic Acid Detection



USAMRIID

Integrated Identification and Diagnostic System




✓ ***Definitive biological agent diagnosis requires integrated identification technologies***



USAMRIID

Biomarkers

Diversity & Depth



Specific virulence markers

Genus and species markers

Common pathogenic markers
& antibiotic resistance

Host Response Markers

✓ **Avoid Technological Surprise**



USAMRIID


Evolutionary Strategy

2002

**Portable Rapid Nucleic
Acid Analysis System**

2003

**Improved
Immunodiagnostic System**



2007

**Comprehensive
Integrated
Diagnostic Systems**




USAMRIID

Medical Diagnostics




Platforms




Assays/Sample Processing



Reagents



Standardization



Validation



USAMRIID

Specimen Processing Challenge

Medical specimens

- swabs
- whole blood and serum
- urine
- feces
- sputum
- lesion exudate
- tissues

Environmental samples


- high volume air sampler
- swabs
- water
- soil

✓ Each matrix may require a unique processing protocol




USAMRIID


Specimen Processing Options




Integrated Cartridge System



Tube and Paper Based Methods



Automated Systems





Manual Cartridge Systems

USAMRIID




Gene Amplification and Detection Option



**Four Cartridge
Integrated Nucleic Acid Analysis System**
USAMRIID



Common Gene Amplification Chemistry



- Common gene amplification chemistry
- Assays for over 26 biological agents
- COTS technology
- Over 50 assays developed



USAMRIID

Gene Amplification Reagents



- Pre-formulated
- Pre-dispensed
- Single-dose
- Ambient-temperature-stable
- Greater reproducibility
- Minimize errors and contamination

Tetracorp, Inc

Dried down
reagents



- Amersham-Pharmacia
Ready-to-Go PCR
PCR-RT Beads
Custom BT Beads

- Idaho Technologies
Lyophilized in
USAMRIID tubes



Leading Instrument Options



Smart Cycler™ XC



RAPID/LightCycler™

- Rugged and portable
- Rapid (25 to 40 mins after specimen processing)
- Sensitive
- Common fluorescent probe chemistry
- Different engineering and operation concepts



USAMRIID

Device Comparisons


Characteristics	ABI 7700 [®]	SmartCycler [®]	R.A.P.I.D.S.
Optical Excitation	1 color	4 color	1 color
Optical Detection	4 color	4 color	2 color
Independent Thermocycling	no	yes	no
Reaction Volume	50 ul	100 ul	20 ul
Rapid Thermocycling	no	yes	yes
Low Power	no	yes	no
Average Analysis Time	2.5 to 3.0 h	20 to 40 min	20 to 40 min



USAMRIID


Evaluation Trials

- Laboratory trials
- Animal models
- Field studies
- Hospital-based



USAMRIID

ECL Immunoassay




- High sensitivity
- Wide dynamic range
- 15 min assay
- Stable reagents



USAMRIID



Electrochemiluminescence (ECL) Reaction First Generation Device - ORIGEN[®]



USAMRIID

ECLM - Igen

- Modular design
- Conducive to portability
- Clinical device uses eight modules
 - 96 well plates
 - ease of maintenance



USAMRIID

Single Tube ECL Assay



USAMRIID


ECL Assays

- SEB Toxin
- Ricin Toxin
- Bot A Toxin
- Plague F1 Antigen
- Anthrax PA
- VEE Virus
- Orthopox
- Near Future
 - Fielding
 - » TAML
 - » Multi-center trials
 - Miniaturization and optimization of next generation device




Magnetic Assay Reader

- Measures the amount of magnetic material in the analytical region
- Oscillating magnetic field used to excite paramagnetic particles
- Detector measures the local magnetic field



USAMRIID

Hand Held Assay Strip



- Simple
- Multiple analyte testing possible
- 10-15 min test




USAMRIID

Comprehensive Integrated Diagnostic System




Automated Sample Processing



Immunodiagnosics




Nucleic Acid Detection



USAMRIID




Far Term Development



Chip Technology



Mass Spectrometry



Advanced Nucleic
Acid Sequencing



USAMRIID

Summary

- Emerging technologies will eventually allow for the development of light weight hand-held devices
- Near term development is focused upon new methods of amplification of gene products and detection of immune complexes
- Far term development may take advantage of microchip, mass spectrometry, and other technologies with broad applications

